




**Accessories and Supports for CON-TREX**




# 5054  
Footplates rigid with snowboard binding (pair)




# 5053  
Footplates flexible with snowboard binding (pair)



# 5050  
Footplates rigid (pair)




# 5049  
Footplates standard (pair)




# 5026  
Strap 45 cm (for 5050, 5049)




# 5014  
Strap 60 cm (for 5050, 5049)




# 5077  
Emergency stop switch



# 5078  
Motion control yellow/blue



# 5009  
Additional dynamometer clamping device



# 5068  
Wall board for standard adapters



All **CON-TREX** Systems conform to the **CE** regulations


Changes in construction, equipment and supply as well as deviation in colour remain subject to alteration

CMV AG, Zürichstrasse 69, CH-8600 Dübendorf, Switzerland


Tel. +41 (0)44-802 12 72, Fax +41 (0)44-802 12 77

Internet: <http://www.con-trex.ch>


e-mail: [info@con-trex.ch](mailto:info@con-trex.ch)




# 5001  
Base adapter long




# 5016  
Base adapter for free movement with 5001 (and 5063 or 5064) and 5060, 5061 or 5062




# 5007  
Knee adapter with 5001




# 5015  
Hip adapter with 5001




# 5055  
ACL-Protection adapter (Anterior cruciate ligament) with 5001




# 5020  
Shoulder adapter Ex/Flex, Add/Abd with 5001



# 5033  
Shoulder rotation with 5001




# 5010  
Support clamping device with 5022 or 5032 and 5011 or 5029




# 5022  
Support rail short with 5010 and 5011 or 5029



# 5032  
Support rail long with 5010 and 5011 or 5029




# 5058  
Stabilizer shoulder rotation



# 5023  
Leg stabiliser for MJ and WS with 5001 and 5011 or 5029



# 5011  
Support narrow with 5010 and 5022 or 5032




# 5029  
Support large with 5010 and 5022 or 5032




# 5063  
Distance rod 30cm for free movement with 5001, 5016 and 5060, 5061 or 5062




# 5064  
Distance rod 60cm for free movement with 5001, 5016 and 5060, 5061 or 5062



# 5060  
Single grip rigid for free movement with 5001 and 5016 (and 5063 or 5064)



# 5061  
Single grip flexible for free movement with 5001 and 5016 (and 5063 or 5064)



# 5062  
Double grip for free movement with 5001 and 5016 (and 5063 or 5064)

# CON-TREX Biomechanical Test- and Training Systems

# CON-TREX Biomechanical Test- and Training Systems

